

## The Toxmatrix: Chemo-genomic profiling identifies interactions that reveal mechanisms of toxicity

Zhi-Bin Tong, Ruili Huang, Yuhong Wang, Carleen Klumpp-Thomas, John Braisted, Zina Itkin, Paul Shinn, Menghang Xia, Anton Simeonov, and David Gerhold

*Chem. Res. Toxicol.*, **Just Accepted Manuscript** • DOI: 10.1021/acs.chemrestox.7b00290 • Publication Date (Web): 20 Nov 2017

Downloaded from <http://pubs.acs.org> on November 25, 2017

### Just Accepted

“Just Accepted” manuscripts have been peer-reviewed and accepted for publication. They are posted online prior to technical editing, formatting for publication and author proofing. The American Chemical Society provides “Just Accepted” as a free service to the research community to expedite the dissemination of scientific material as soon as possible after acceptance. “Just Accepted” manuscripts appear in full in PDF format accompanied by an HTML abstract. “Just Accepted” manuscripts have been fully peer reviewed, but should not be considered the official version of record. They are accessible to all readers and citable by the Digital Object Identifier (DOI®). “Just Accepted” is an optional service offered to authors. Therefore, the “Just Accepted” Web site may not include all articles that will be published in the journal. After a manuscript is technically edited and formatted, it will be removed from the “Just Accepted” Web site and published as an ASAP article. Note that technical editing may introduce minor changes to the manuscript text and/or graphics which could affect content, and all legal disclaimers and ethical guidelines that apply to the journal pertain. ACS cannot be held responsible for errors or consequences arising from the use of information contained in these “Just Accepted” manuscripts.

# The Toxmatrix: Chemo-genomic profiling identifies interactions that reveal mechanisms of toxicity

Zhi-Bin Tong, Ruili Huang, Yuhong Wang, Carleen A. Klumpp-Thomas, John C. Braisted, Zina Itkin, Paul Shinn, Menghang Xia, Anton Simeonov and David L. Gerhold\*


National Center for Advancing Translational Sciences (NCATS), National Institutes of Health, 9800 Medical Center Drive C-345E, Bethesda, Maryland 20892

\*Corresponding Author

Email address: [david.gerhold@nih.gov](mailto:david.gerhold@nih.gov)

phone: 301-827-5346

address: 6701 Democracy Boulevard, Bethesda MD 20892-4874



## ABSTRACT

A chemical genomics ‘Toxmatrix’ method was developed to elucidate mechanisms of cytotoxicity using neuronal models. Quantitative high-throughput screening (qHTS) was applied to systematically screen each toxicant against a panel of 70 modulators, drugs or chemicals that act on a known target, to identify interactions that either protect or sensitize cells to each toxicant. Thirty-two toxicants were tested at 10 concentrations for cytotoxicity to SH-SY5Y human neuroblastoma cells, with results fitted to the Hill equation to determine an  $IC_{50}$  for each toxicant. Thirty-three toxicant:modulator interactions were identified in SH-SY5Y cells for 14 toxicants, as modulators that shifted toxicant  $IC_{50}$  values lower or higher. The target of each modulator that sensitizes cells or protects cells from a toxicant suggests a mode of toxicant action or cellular adaptation. In secondary screening, we tested modulator-toxicant pairs identified from the SH-SY5Y primary screening for interactions in three differentiated neuronal human cell lines: dSH-SY5Y, conditionally-immortalized dopaminergic neurons (LUHMES), and Neural Stem Cells. Twenty toxicant-modulator pairs showed pronounced interactions in one or several differentiated cell models. Additional testing confirmed that several modulators acted through their primary targets. For example, several chelators protected differentiated LUHMES neurons from four toxicants by chelation of divalent cations; and buthionine sulphoximine sensitized cells to 6-hydroxydopamine and 4-(Methylamino)phenol hemisulfate by blocking glutathione synthesis. Such modulators that interact with multiple neurotoxicants suggest these may be vulnerable toxicity pathways in neurons. Thus, the Toxmatrix method is a systematic high-throughput approach that can identify mechanisms of toxicity and of cellular adaptation.

**Keywords:** neurodegeneration, high-throughput screening, neuronal models

## INTRODUCTION

Toxic environmental chemicals that have traditionally been identified by animal testing, are increasingly being identified by *in vitro* assays. For example, *in vitro* quantitative high throughput screening (qHTS) has been performed by the Tox21 Consortium for ten thousand such chemicals.<sup>1-3</sup> Specialized qHTS methods were developed by the Tox21 partners to robotically screen thousands of compounds at multiple concentrations in 1536-well microplates to identify those that activate a variety of cellular-stress events. Chemicals that activate such events, *e.g.* disruption of mitochondrial membrane potential, TP53 activation, cytotoxicity, *et al.*,<sup>4</sup> demand further evaluation as putative toxicants. We evaluated whether one could extend this qHTS approach to systematically elucidate mechanisms for toxicants identified by qHTS by identifying interactions between each toxicant and toxicity pathways within the cell model. Thus, toxicants were tested using a systems-biology approach for interaction with a panel of modulators, chemicals that each perturb a cell stress-related pathway, to determine if they protect or sensitize cells to each toxicant.

Many drugs and drug metabolites are selectively toxic to a specific human cell type, either through selective pharmacokinetic properties (PK), or selective toxicodynamics (TD).<sup>5</sup> PK may determine selective toxicity to neurons, *e.g.* when MPTP is metabolized into toxic MPP<sup>+</sup> by astrocytes.<sup>6</sup> In turn, the “Toxicity Pathways” approach<sup>7</sup> describes TD toxicity mechanisms as normal cellular pathways that result in adverse events when disrupted. This approach implies that the sensitive cell type has a unique pathway, or is especially sensitive to disruption of the toxicity pathway. Smirnova *et al.* have suggested that cellular resilience, the manner in which a cell adapts to resist toxicity, may be more important than the cells’ direct susceptibility to toxicity in determining which cell types survive.<sup>8</sup> The toxicity pathway and resilience approaches are not conflicting, but rather complementary descriptions of how the toxicant insults the cell, and how the cell attempts to adapt to the insult, respectively.

Neurons are protected from most toxicants by the blood-brain barrier, yet neurons must have unique TD pathways that make these cells particularly sensitive to toxicants. For example, neurotoxicants may be

1  
2  
3 toxicodynamically selective for neurons by disrupting the microtubules that support the long axons of  
4 neurons, or by disrupting mitochondria that are used intensively to repolarize membranes in electrically-  
5 active neurons. Characterization of neurotoxicants that kill mature neurons is urgently needed in order to  
6 identify neurotoxic chemicals that are suspected to cause Parkinson's disease (PD), Alzheimer's disease  
7 (AD), and other neurodegenerative disorders.<sup>9</sup>

8  
9  
10  
11  
12  
13  
14 Neuronal cytotoxicity models were investigated in a previous study from this laboratory that compared 3  
15 human neuronal models: SH-SY5Y neuroblastoma cells,<sup>10</sup> LUHMES conditionally-immortalized  
16 dopaminergic neurons,<sup>11</sup> and Neural Stem Cells (NSC)<sup>12</sup> that differentiate into a mix of neuronal,  
17 astrocytic, and oligodendrocytic cell types.<sup>13</sup> For neurotoxicological screening the LUHMES model was  
18 preferred due to its ability to differentiate nearly 100% of cells within seven days, to its high-level  
19 expression of neuronal markers, and to its greater sensitivity to 21 -out-of- 32 known or candidate  
20 neurotoxicants tested.<sup>12</sup> These three models, and particularly dLUHMES cells, are further applied in this  
21 study to elucidate the mode or mechanisms by which neurotoxicants may kill neurons. Thus, LUHMES  
22 cells seem to recapitulate *in vitro* the selective sensitivity of *in vivo* neurons to neurotoxicants in many  
23 cases. Nevertheless, some neurotoxicants will demand a model system that includes metabolic activation,  
24 neuronal interactions with other cell types, or the 3-dimensional microenvironment found within the  
25 brain.

26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41 Toxicants have previously been characterized to elucidate their mechanisms of action, or classify modes  
42 of toxicity, by identifying a second chemical factor that either sensitizes or protects the cell from death.  
43  
44  
45 Wolpaw *et al.* demonstrated a systematic "modulatory profiling" screen to identify pairs of drugs that  
46 showed synergy or antagonism in killing tumor-derived cells.<sup>14</sup> Drugs that showed similar modulatory  
47 profiles were grouped in order to find outlier anticancer drug candidates that implied a novel mode of  
48 action.

1  
2  
3 In the present study we evaluate a systematic and quantitative matrix screening method called the  
4 Toxmatrix, to determine whether it can yield insights into toxicant mechanisms. The Toxmatrix method  
5 uses modulator chemicals that address known targets, to systematically identify those that promote or  
6 inhibit toxicity of known or suspected neurotoxicants. Seventy modulators were selected from the  
7 published literature to target pathways for cellular adaptation to stress, or for pathways that regulate  
8 necrosis or apoptosis, cell survival, or regulation of energy metabolism. In brief, master plates for  
9 toxicants and modulators are prepared using robots and barcoding to enable automated sample-tracking.  
10 qHTS uses cells cultured in 5  $\mu$ L volumes in 1536-well microplates, dosed with precise 100 nL volumes  
11 using an acoustic dispenser. Concentration-response cytotoxicity data for each toxicant and toxicant +  
12 modulators are recorded in a database and automatically analyzed by curve-fitting, derivation of IC<sub>50</sub>  
13 values and quality scores, and visualization. This qHTS infrastructure enabled us to formulate and  
14 systematize the Toxmatrix method and to evaluate its suitability to determine mechanisms of cytotoxicity.  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28

## 29 **METHODS**

30  
31  
32 **Chemicals and Compound Management** The 32 toxicant chemicals and 70 modulators were  
33 included in this study were purchased from Sigma-Aldrich Chemical Co., Selleck Chemicals LLC, R &  
34 D/Tocris, or synthesized at NCATS (see supplementary Table 1). The compounds were dissolved in  
35 dimethyl sulfoxide (DMSO) at either 20 mM (200X) or 50 mM (500X) for the toxicants. Modulators  
36 were diluted in 3 five-fold steps to generate 200x master solutions for primary screening in SH-SY5Y  
37 cells, to include a middle dosage selected from published studies. For secondary screening in three cell  
38 lines, toxicants were diluted to 9 concentrations in three-fold steps and modulators were also diluted to 9  
39 concentrations in three-fold steps based on the minimum toxic concentration from primary screening.  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49

50 **Cell Culture and Differentiation** Human cell lines including neuroblastoma cell line SH-SY5Y,  
51 neuronal stem cell NSC, and a conditionally-immortalized human fetal mesencephalic cell line LUHMES  
52 cells were cultured and differentiated as described before.<sup>6, 12</sup> Briefly, SH-SY5Y cells (ATCC,  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 Gaithersburg, MD) were cultured in 1:1 mixture (EME/F12) of Eagle's Minimum Essential (EMEM,  
4 ATCC) and F12 (Invitrogen) containing 10% fetal bovine sera (FBS). SH-SY5Y cell differentiation was  
5 induced by 10  $\mu$ M all-trans-retinoic acid in 3% FBS culture medium. Gibco® Human Neural Stem Cells  
6 (NSC; Invitrogen/Thermo-Fisher, Waltham, Massachusetts, USA) were adherently cultured in vessels  
7 coated with Geltrex and grown in the defined StemPro NSC SFM complete medium consisting of  
8 KnockOut D-MED/F-12 with 2% StemPro Neural Supplement, 2 mM GlutaMax-I Supplement, and  
9 human recombinant bFGF and EGF (20 ng/ml each). NSC differentiation was carried out in StemPro  
10 NSC SFM medium without bFGF and EGF.<sup>12</sup> LUHMES cells were grown in culture vessels coated with  
11 poly-L-ornithine and human fibronectin (50  $\mu$ g/ml each) in Advanced DMEM/F12 medium containing 2  
12 mM L-glutamine, N2-supplement (1X) and 40 ng/ml bFGF. LUHMES cells were subsequently  
13 differentiated by adding 1  $\mu$ M cAMP and 1  $\mu$ g/ml tetracycline in advanced DMEM/F12 medium with 2  
14 mM L-glutamine, 1X N2-supplement, and 2ng/ml GDNF.<sup>6</sup> Both proliferating and differentiating cells  
15 were transferred into tissue culture-treated 1536-well plates with the same coatings using a Multidrop  
16 instrument (Thermo-Fisher Scientific, Waltham, MA USA) in a volume 5  $\mu$ L containing 2000 cells per  
17 well at one day before the modulators and toxicants were applied. After 9 passages, cell lines were  
18 discarded and replaced from frozen stocks. Identities of the three cell lines were verified following each  
19 experiment by short tandem repeat profiling (WiCell Research Inst. Madison, Wisconsin, USA). Cell  
20 lines were also checked following each experiment for mycoplasma contamination using the  
21 MycoAlert™ Mycoplasma Detection Kit (Lonza, Basel Switzerland).

22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46 **qHTS Toxmatrix Screening** Undifferentiated cells, or cells allowed to differentiate in T225 flasks,  
47 were distributed in 1536-well plates one day before applying chemical treatments, using 2,000 cells per  
48 well in 5  $\mu$ L medium. The primary Toxmatrix consisted of 32 chemical toxicants and 70 modulators  
49 (supplementary Table 1). The toxicants and modulators were made in separate 1536-well compound  
50 plates in appropriate concentrations as 200X master plates. Each of the toxicants was applied to cells at 10  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 concentrations and no-toxicant controls from a maximum of 100  $\mu\text{M}$  in three-fold steps; while each of the  
4 modulators was applied to cells at 3 concentrations in five-fold steps, designated “low”, “middle” and  
5 “high”. For this “10x3” primary screening, each compound was loaded in 25 nL of DMSO using a 1536  
6 pintool instrument (V&P Scientific). In the secondary 9x9 Toxmatrix screening, 36 pairs of toxicants (13)  
7 and modulators (22) with appropriate controls were selected, based on data from the primary screening.  
8 Each compound was serially diluted (1:3) to make compound source plates at 500-fold the final  
9 concentration, including 9 concentrations of the toxicants and 9 concentrations for the modulators, plus no  
10 toxicant and no-modulator samples. Each dosage was loaded in 10 nL into cultured cells using an ATS-  
11 100 acoustic dispenser (EDC Biosystems, Fremont, CA, USA). In both primary and secondary screenings  
12 cells were cultured in 5  $\mu\text{L}$  of medium in 1536-well plates. The modulators were loaded four h before the  
13 toxicants were applied, and every permutation of modulator- and toxicant concentrations was tested. The  
14 first 4 rows of cells in each 1536-well plate were treated with DMSO only as negative controls for cell  
15 survival, or tetra-octyl ammonium bromide as positive controls for 100% cell killing. Cytotoxicity assays  
16 were carried out at 24 and 48 h after the modulators and toxicants were applied.  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32

33  
34 **Cell Viability and Caspase 3/7 Assays** Cells were assessed for viability using a CellTitreGlo®  
35 luminescence assay, Promega, Madison Wisconsin, USA). After compound treatment, 5  $\mu\text{L}$  of these  
36 detection reagents were added into each well of the 1536-well assay plates, and the assay plates were  
37 incubated for 30-60 min at room temperature. Luminescent intensity was measured on a Viewlux  
38 instrument (PerkinElmer, Waltham MA, USA) to indicate cell viability.  
39  
40  
41  
42  
43  
44

45 **qHTS Data Analysis** Analysis of compound concentration–response data was performed as  
46 previously described.<sup>4, 15</sup> Briefly, raw plate reads were first normalized relative to the positive control  
47 compound (-100%) and DMSO-only wells (0%) as follows: % Activity =  $((V_{\text{compound}} - V_{\text{DMSO}})/(V_{\text{pos}} -$   
48  $V_{\text{DMSO}})) \times 100$ , where  $V_{\text{compound}}$  denotes the compound well values,  $V_{\text{pos}}$  denotes the median value of the  
49 positive control wells, and  $V_{\text{DMSO}}$  denotes the median values of the DMSO-only wells, and then corrected  
50 by applying a NCATS in-house correction algorithm.<sup>16</sup> Concentration–response titration points for each  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 toxicant at fixed modulator concentration were fitted to a four-parameter Hill equation<sup>17</sup> using an iterative  
4 grid algorithm.<sup>18</sup> These curve-fitting results yielded concentrations of half-maximal inhibition (IC<sub>50</sub>) and  
5 maximal response (efficacy) values. Modulator concentrations that caused  $\geq 10\%$  cytotoxicity were  
6 excluded. Modulator:toxicant:cell-type combinations that yielded a significant IC<sub>50</sub> shift, exceeding 3-fold  
7 either leftward (+), or rightward (-) were verified by visual inspection from the curve fits. The maximum  
8 IC<sub>50</sub> shift is reported among non-cytotoxic modulator concentrations.  
9  
10  
11  
12  
13  
14  
15

16 **Validation of Toxmatrix Screenings** Further experiments to validate the selected pairs of toxicant-  
17 modulator functional interactions were carried out manually in 384 well plates using differentiated human  
18 LUHMES cells. Briefly, 10,000 LUHMES cells/well were loaded in 50  $\mu\text{L}$  to 384 well plates and  
19 differentiation was induced for 7 days. Modulators PHE, M30, or BSO were used to treat cells for a four h  
20 interval prior to treatment with a toxicant concentration that caused approximately 50% killing.  
21 Supplements GSH or NAC were also added for this 4 hour interval as indicated. After treatment, cell  
22 viability was assayed at 24 and 48 h to determine whether a modulator or modulator + supplement  
23 affected toxicant-induced cytotoxicity.  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33

## 34 35 36 37 **RESULTS AND DISCUSSION**

38 **Identification of Toxicant:Modulator Pairs in SH-SY5Y Cells** We sought to build a matrix  
39 screening tool that can reveal mechanisms of toxicity or adaptation by identifying modulators that protect  
40 or sensitize cells to the toxicant. Seventy modulators were selected from the literature to target pathways  
41 related to necrosis or apoptosis, cell survival, regulation of energy metabolism, as well as adaptation to  
42 stress including: detoxification, oxidative stress, DNA damage and repair, cell cycle regulation, innate  
43 immune responses, ER stress, autophagy, proteolysis, and disruption of enzyme cofactors. This  
44 Toxmatrix approach was designed to identify modulators for each toxicant that shifted the concentration-  
45 response curve of a toxicant left or right relative to the modulator-free controls (sensitization or  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 protection, respectively). The goal was to identify a modulator concentration that shifted the  $IC_{50}$   
4 maximally without exhibiting cytotoxicity by itself. Cells were allowed to respond to modulator  
5 chemicals for four hours prior to toxicant treatments. SH-SY5Y cells were used initially to screen of 32  
6 toxicants x 70 modulators. The SH-SY5Y cell line neuronal model was selected because it was used to  
7 identify several candidate neurotoxicants,<sup>1, 4, 12</sup> and because SH-SY5Y cells were readily cultured for  
8 high-throughput screening<sup>1</sup> in 1536-well microplates. The qHTS screening format of Inglese *et al.*<sup>19</sup> was  
9 adopted to assay cytotoxicity in 1536-well microplates using 10 concentrations of each toxicant from 100  
10  $\mu$ M to 5 nM in three-fold increments, and no-toxicant controls; as well as 70 modulators at three  
11 concentrations each, and no-modulator controls, for each query toxicant.

12  
13 This survey using SH-SY5Y cells examined combinations of 32 toxicants at 10 concentrations each (plus  
14 the vehicle no-toxicant control), with 70 modulators at 3 concentrations each (plus modulator-only  
15 controls); comprising 98,560 nominal sample-data points. For each toxicant:modulator pair, a curve was  
16 drawn to fit the Hill equation<sup>18</sup> for the 10 toxicant concentrations without modulator, and an additional  
17 curve for each modulator concentration. Representative plots are shown for three toxicant:modulator  
18 pairs effects on SH-SY5Y cells in Figure 1. Figure 1 shows the cytotoxic concentration-response for the  
19 toxicant MAP, left-shifted hence sensitized by modulator  $CoCl_2$  in Fig1A, unchanged by treatment with  
20 modulator necrostatin 5 in Fig 1B, and right-shifted hence protected by treatment with modulator NAC in  
21 Fig 1C. Among 64 microplates, the median % Coefficient of Variation for DMSO controls was 3.9% with  
22 a standard deviation (SD) of 0.71%, the median signal/background ratio was 31 with a SD of 1.7, and the  
23 median  $Z'$  score<sup>4</sup> was 0.90 with a SD of .02. These reproducible data and robust curve-fits resulted in  
24 reliable automated classification of protection or sensitization events for nearly all toxicant:modulator  
25 pairs. Modulators that caused cytotoxicity without a toxicant at high concentrations (viability <90%, *e.g.*  
26 the high concentration of  $CoCl_2$  in Fig 1A) were excluded from results. In several cases automatic calls of  
27 sensitization or protection appeared to result from an upward- or downward- shift exceeding 10% relative  
28 to the 100% baseline for particular modulators at low toxicant concentration. Such curves were also

1  
2  
3 excluded from results as likely artefacts. This process enabled us to identify discrete and substantial  
4 protective or sensitizing interactions for further studies, and to exclude false-positive or weak interactions.  
5  
6

7  
8 Figure 2 summarizes 33 significant interactions among pairwise combinations of 32 toxicants with 70  
9 modulators, in which the modulator increased toxicity (sensitization) in 23 pairs, and the modulator  
10 decreased toxicity (protection) in 10 pairs. The significant toxicant:modulator pairs indicated in Figure 2  
11 were evaluated by examining the curves for consistent sensitization or protection by the modulator pairs  
12 (modulator caused >3-fold shift of the IC<sub>50</sub> at a concentration that caused <10% cytotoxicity) as  
13 illustrated in Fig 1A and Fig 1B. Figure 2 indicates several generalizable results. First, a particular  
14 modulator tended to function either in protection or sensitization, but rarely in both modes; whereas  
15 several toxicants were subject to both effects by modulators. Among the 70 modulators tested, five were  
16 protective whereas sixteen acted as sensitizers. Two modulators functioned in both modes, the divalent  
17 cation chelator DMPS, which protected SH-SY5Y cells from 5 toxicants, but sensitized cells to HCP; and  
18 KU0063794, a TOR inhibitor drug that both protected SH-SY5Y cells from OTZ, and sensitized them to  
19 MHG. Unexpectedly, several modulators that were designed to protect cells from TNFalpha-driven  
20 necrosis, necrostatins 1, 5, and 7, instead sensitized SH-SY5Y cells to killing by HCP or MALG under  
21 these conditions (Figure 2).  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37

38 **Toxicant:Modulator Interactions in 3 Neuronal Models** From a matrix of 32 toxicants and 70  
39 modulators, we identified 33 pairwise interactions in SH-SY5Y cells (Figure 2), including combinations  
40 of 14 known or suspected neurotoxicants and 19 modulators. These 33 pairs were selected for secondary  
41 screening using SH-SY5Y, LUHMES, and NSC cells, after differentiation protocols that increased the  
42 neuronal characteristics of these models.<sup>12</sup> Conditions were developed for each cell line, extending  
43 previous work using 384-well microplates,<sup>12</sup> to enable differentiation and culture in 1536-well format for  
44 high-throughput Toxmatrix cytotoxicity assays. In these experiments, cells were differentiated for 7  
45 days,<sup>12</sup> plated, and allowed to re-elaborate neurites (termed dSH-SY5Y, dLUHMES, and dNSC), before  
46 addition of modulators and toxicants. These were performed in 1536-well microplates using 9  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 concentrations of each toxicant and no-toxicant DMSO controls, as well as 9 concentrations of each  
4 modulator and no-modulator controls. For Figure 3, cytotoxicity assays were performed at both 24 h and  
5 48 h to observe both transient interactions that manifested at 24 h and slowly developing interactions that  
6 manifested at 48 h. Cytotoxicity results from this matrix of 9,240 data points were then examined to  
7 identify combinations of toxicants and modulators that interacted selectively and robustly in these  
8 differentiated cells. The maximum fold-shift in  $IC_{50}$  is reported in Figure 3 for each interaction that  
9 resulted in a significant shift in  $IC_{50}$  (>3-fold increase or decrease) at a modulator concentration that  
10 resulted in <10% cytotoxicity measured as intracellular [ATP].  
11

12  
13 Figure 3 illustrates two notable trends in toxicant:modulator:cell line interactions. Each modulator  
14 characteristically either protected or sensitized cells with multiple toxicants; and toxicant:modulator  
15 interactions tended to occur in more than one model. Specifically DMPS, NAC, nicotinamide, and SP-  
16 600125 showed only protective actions, whereas  $CoCl_2$ , BSO, necrostatin 7, WYE-354, CNF-2024,  
17 SRT1720, GSK-1904529A and TW-37 showed only sensitizing actions. DMPS protected cells from five  
18 toxicants, MAP, 6HD, CAP, DTCM, and MHG; and N-acetyl-L-cysteine (NAC) protected cells from  
19 MAP, DTCM, and MHG. Conversely, necrostatin 7 sensitized cells to three toxicants: 6HD, MALG, and  
20 MHG. In fact, several modulators exhibited consistent effects for multiple toxicants, particularly DMPS,  
21 NAC, and necrostatin 7. Similarly, interactions almost invariably occurred in treatments with the same  
22 toxicant in multiple cell types, the exceptions being GSK-1904529A that sensitized only dSH-SY5Y from  
23 MALG, and nicotinamide that protected only dLUHMES from VIN. Note that in Figure 3  $IC_{50}$  fold  
24 change is indicated for changes exceeding 3-fold, whereas changes between 2-3-fold are colored red or  
25 blue without a number to reveal trends that did not reach the 3-fold significance standard. Interactions  
26 were often but not always observed at both 24 h and 48 h time points. For example, at both time points  
27 DMPS protected: dLUHMES from 6HD, all three cell lines from CAP, and dSH-SY5Y as well as  
28 dLUHMES from DTCM. Conversely, enhanced cytotoxicity was sometimes observed at only one time  
29 point; *e.g.* dLUHMES was selectively sensitized to 6HD by both BSO and necrostatin 7 at 24 h but only  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 weakly or not at all at 48 h. Many toxicant:modulator:cell line interactions gave rise to testable  
4  
5 hypotheses; these are discussed below.  
6  
7

8  
9 **Sensitizing Interactions in Differentiated Cell Models** Toxicant:modulator:cell line interactions that  
10 sensitized neuronal cells implied that the modulator acted upon a target-pathway that enhanced toxicity or  
11 interfered with adaptation. The diversity of significant toxicant:modulator combinations generated a  
12 variety of hypotheses regarding toxicant modes of action. For example, sensitization of dLUHMES and  
13 dNSC to MAP by  $\text{CoCl}_2$ , an inducer of the HIF1 hypoxia pathway, suggests that HIF1 activation  
14 stimulated a pro-apoptotic or pro-necrotic phenotype that made dLUHMES susceptible to MAP cytotoxic  
15 challenge. Necrostatin 7 sensitized several cell lines to 6HD, MALG, and MHG. This was surprising,  
16 since the necrostatins were developed for their ability to protect various cell types from the cytotoxic  
17 effect of Tumor Necrosis Factor.<sup>20</sup> There is precedent for the presumed target of the necrostatins, RIP1K,  
18 acting on the NF $\kappa$ B pathway, either sensitizing or protecting cells dependent on cell type and cytotoxic  
19 stimuli.<sup>21</sup> Nevertheless, it was shown recently that necrostatin 7 does not inhibit RIP1 Kinase, the known  
20 target of necrostatin 1, and -5, hence necrostatin 7 may work through a different, unknown, target.<sup>20</sup> In  
21 support of the possibility that necrostatin 7 works *via* a unique target, necrostatins 1 and 5 coincided in  
22 sensitizing SH-SY5Y cells to HCP; whereas necrostatin-7 differed in sensitizing this cell line to MALG  
23 and CR2. The sensitizing activity of necrostatin 7 to 6HD, MALG, and MHG in these neuronal models  
24 suggested a shared pathway of toxicity among these three toxicants, but the mode or mechanism awaits  
25 further characterization. BIIB021, an inhibitor of HSP90 that inactivates NF $\kappa$ B,<sup>22</sup> also sensitized dSH-  
26 SY5Y and dLUHMES to CLM. SRT1720, a Sirtuin activator, sensitized dSH-SY5Y, dLUHMES, and  
27 dNSC cell lines to MPP. This was surprising given that SIRT1 activators as well as SIRT2 inhibitors have  
28 been proposed as therapeutic treatments for PD.<sup>23</sup> Since SRT1720 has been shown to increase  
29 mitochondrial metabolism, and MPP is thought to attack mitochondrial electron transport, we hypothesize  
30 that SRT1720 increases electron flow into an electron transport chain that is blocked at Complex I by  
31 MPP, resulting in increased electron flow to  $\text{H}_2\text{O}_2$  and increased toxic radical formation.<sup>24</sup> GSK-  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 1904529A, an IGF1 Receptor inhibitor, was observed to sensitize dSH-SY5Y cells to MALG, an  
4 antifungal used in aquaculture; possibly by blocking unknown survival factor(s) that are dependent on  
5 IGF1 action.<sup>25</sup> TW-37, an inhibitor of the antiapoptotic BCL2 protein, sensitized dSH-SY5Y cells to  
6 TUNC, suggesting that BCL2 mitigates an apoptotic activity of TUNC.<sup>26</sup>  
7  
8  
9

10  
11  
12 **Protective Interactions in Differentiated Cell Models** A variety of modulators were identified that  
13 protected cells from known or suspected neurotoxicants, including NAC and DMPS. Because NAC and  
14 DMPS protected neuronal models from a variety of neurodegenerative toxicants, interactions involving  
15 these two modulators were selected for further studies aimed at confirming or refuting the expected  
16 targets of NAC and DMPS. NAC protected dSH-SY5Y from MAP, DTCM, and MHG, as well as  
17 protecting dLUHMES from DTCM (Figure 3). This protective activity indicates that these toxicants cause  
18 oxidative stress, originating either from reactive drug metabolites or oxidative byproducts of H<sub>2</sub>O<sub>2</sub>  
19 resulting from disrupted electron transport in the mitochondria.<sup>27</sup> NAC can act either by feeding GSH  
20 synthesis or by acting directly as a scavenger, a question that is addressed below in Figure 4. DMPS also  
21 protected cells from multiple toxicants; protecting dSH-SY5Y and dNSC from MAP, dLUHMES and  
22 dNSC from 6HD and MHG, dSH-SY5Y and dLUHMES from DTCM, and all three cell lines from CAP,  
23 Figure 3. These interactions reveal that cytotoxic mechanisms of OTZ, MAP, 6HD, CAP, DTCM, and  
24 MHG involve liberation of divalent cations, likely to be free iron, which is known to generate toxic  
25 hydroxyl radical from H<sub>2</sub>O<sub>2</sub>.<sup>28</sup> Nicotinamide protected dLUHMES cells from the toxicity of MHG, acting  
26 either as a SIRT inhibitor, or as an antioxidant. JNK 1/2/3 inhibitor SP-600125 mitigated the toxicity of  
27 VIN towards both dLUHMES and dNSC cells, suggesting that the JNK stress pathway functions as a  
28 toxicity pathway downstream of VIN treatment.  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48

49 Modulator interactions with toxicants 6HD, CAP, MAP, and DTCM were of particular interest to study in  
50 dLUHMES cells, since dLUHMES are more sensitive to these toxicants than either undifferentiated  
51 LUHMES, or SH-SY5Y cells.<sup>12</sup> These results suggested that toxicants 6HD, CAP, MAP, and DTCM may  
52 be selectively neurotoxic because of neuronal- or dopaminergic-neuronal characteristics of the  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 dLUHMES cells. dSH-SY5Y and dLUHMES cells were protected from these toxicants by DMPS (except  
4 6HD in SH-SY5Y and MAP in dLUHMES did not meet the 3-fold threshold), which led us to evaluate  
5 DMPS action using additional chelators. Additionally, dLUHMES and dSH-SY5Y were protected from  
6 DTCM by NAC and sensitized to 6HD by BSO, which causes depletion of glutathione, which led us to  
7 additional study of the role of glutathione in adaptation to oxidative stress. These additional evaluations of  
8 BSO and DTCM were performed in dLUHMES neurons to determine whether interactions arose from the  
9 generally-accepted target of the specific modulator, revealing pathways impacted by each toxicant, or  
10 arose from off-target activities that obscure these pathways. In these experiments, fixed concentrations of  
11 toxicant and modulator were selected to cause a level of cytotoxicity intermediate between 0% and 100%,  
12 in order to look for increased or decreased toxicity in response to either a supplement or a substitute  
13 modulator.  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25  
26  
27

28 **The Glutathione Pathway in dLUHMES** BSO markedly sensitized dLUHMES cells, as well as  
29 dNSC and dSH-SY5Y, to 6HD (Figure 3). Since BSO is known to block glutathione biosynthesis by  
30 inhibiting  $\gamma$ -glutamylcysteine synthetase,<sup>29</sup> we hypothesized that cells were using GSH to detoxify 6HD,  
31 and that BSO interfered with this detoxification by depleting cellular glutathione. We tested this  
32 hypothesis by a simple supplementation strategy. dLUHMES cells were supplemented with GSH or with  
33 the cysteine precursor NAC, to ask whether these supplements restore the dLUHMES cells' ability to  
34 resist cytotoxicity from 6HD. Figure 4A shows that NAC or GSH partially alleviated the toxicity of 6HD  
35 + BSO to dLUHMES cells, supporting the hypothesis that GSH allays the cytotoxicity of 6HD. GSH at  
36 100  $\mu$ M was particularly effective at decreasing cytotoxicity at an intermediate BSO concentration of 5  
37  $\mu$ M. In cells, 6HD is thought to give rise to a quinone, with production of  $H_2O_2$ . Previous work has  
38 demonstrated that 6HD applied to SH-SY5Y cells is converted into a quinone that gives rise to  $H_2O_2$ ,  
39 causing oxidative stress.<sup>30</sup> The quinone is detoxified by conjugation to glutathione<sup>27, 30</sup> such that GSH is  
40 consumed and regenerated in this process.<sup>35, 38, 31</sup> Our result is also in agreement with a previous  
41 experiment that demonstrated GSH can partially mitigate toxicity of 6HD to SK-N-SH neural cells, the  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 precursor to SH-SY5Y cell.<sup>27</sup> Similarly, BSO increased the sensitivity of dLUHMES cells to MAP, and  
4 supplementation with NAC or GSH partly alleviated this sensitivity (Figure 4B). In the cases of 6HD and  
5 MAP, NAC likely acted at least partly by replenishing GSH, since GSH supplementation was  
6 demonstrated to partly restore the sensitizing activity of BSO (Figure 4). Since GSH is not known to  
7 penetrate cells intact, it is not clear whether GSH acted inside cells, or exerted antioxidant activity  
8 extracellularly. Taken together, these results suggest that GSH mitigates cytotoxicity by relieving  
9 oxidative stress caused by both 6HD and MAP.  
10  
11  
12  
13  
14  
15  
16  
17

18 **Protective Divalent Cation Chelation in dLUHMES** The action of DMPS, a relatively cell-  
19 permeable chelator of divalent cations,<sup>32</sup> in protecting dLUHMES cells from 6HD, MAP, CAP, DTCM,  
20 and MHG, was also further evaluated. Patients diagnosed with heavy metal poisoning are often treated  
21 with a chelator such as DMPS to mobilize accumulated toxic metals; however it is possible that DMPS  
22 also serves as a scavenger of reactive oxidants.<sup>33</sup> To resolve this question additional divalent cation  
23 chelators, 1,10,-phenanthroline (PHE) and M30<sup>34</sup> were tested to determine whether they can substitute for  
24 DMPS to protect dLUHMES cells from two of these five toxicants. Interestingly, chelator PHE, like  
25 DMPS, significantly mitigated the toxicities of 6HD and of MAP to dLUHMES cells (Figure 5A and  
26 Figure 5B, respectively); whereas chelator M30 increased the toxic effects of both 6HD and MAP. Since  
27 M30 is designed to be a scavenger of radicals as well as a chelator, but did not protect LUHMES cells,  
28 these results suggest that it is the metal chelating activity of DMPS that protects LUHMES cells from  
29 toxicants, rather than any oxidant scavenger activity. Although all three chelators can tie up divalent  
30 cations such as  $\text{Fe}^{2+}$  or  $\text{Zn}^{2+}$ , such chelators differ in their abilities to penetrate cell membranes and  
31 partition into subcellular organelles as well as their affinities for ligands, hence may have different effects  
32 in intoxicated cells. Interestingly, whereas DMPS also protected dLUHMES from these five compounds  
33 in these follow-up experiments, neither PHE nor M30 protected dLUHMES cells from 20  $\mu\text{M}$  CAP, 4  $\mu\text{M}$   
34 DTCM, or 1  $\mu\text{M}$  MHG (not shown). Previous work has shown that chelators can protect neuronal models  
35 from 6HD; however different chelators were protective in different cell models. For example, Kobayashi  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 *et al.* showed that the widely studied chelator desferroxamine was protective for SH-SY5Y neurons  
4 treated with 6HD.<sup>35</sup> Desferroxamine and M30 have also shown neuron-protective ability in 6HD-treated  
5 rat models;<sup>34, 36</sup> however neither desferroxamine nor M30 were protective in the dLUHMES model (M30  
6 Figure 5A, desferroxamine not shown). Chelators have been tested in Parkinson's disease patients with  
7 intriguing, but not definitive results, and relatively narrow therapeutic margins.<sup>37, 38</sup> Thus, chelators are  
8 diverse in their actions on different toxicants, revealing in this study that the five toxicants 6HD, MAP,  
9 CAP, DTCM, and MHG are related but not identical in their cytotoxic mechanisms.  
10

11 It is also useful to note that several modulators that yielded interactions in primary screening had “sister  
12 modulators”, that is, modulators with the same primary target (Supplementary table S2). For example,  
13 necrostatin 1 and necrostatin 5 both inhibit RIP1K, and both sensitized SH-SY5Y cells to HCP; and  
14 DMPS and PHE both protected dLUHMES from 6HD and MAP. Other sister modulators failed to yield  
15 similar interactions, *e.g.* mTOR inhibitors GSK-1904529A, KU0063794, and WYE-354; HSP90  
16 inhibitors BIIB021 and Alvespimycin; or Nitric Oxide Synthase inhibitors L-NAME and L-NMMA  
17 (Figure 2 and Table S2); suggesting that modulators are often promiscuous or divergent in their  
18 pharmacological properties.  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35

36 **Toxicants with Shared Mechanisms** Toxicants may be grouped based on being sensitized or  
37 protected by shared modulators, suggesting shared mechanisms of toxicity. For example, DMPS protected  
38 dLUHMES cells from 6HD, MAP, CAP, DTCM, and MHG (Figure 3), indicating a shared mechanism of  
39 toxicant action *via* release of divalent cations; likely iron(II).<sup>28</sup> Similarly, BSO sensitization of dLUHMES  
40 cells to 6HD, and MAP (Figure 4) revealed that these toxicants deplete reduced glutathione that can be  
41 partly ameliorated by NAC (which also protected dLUHMES from DTCM, Figure 3) or by GSH  
42 supplements. Taken together, these results suggest similar modes of toxicity for 6HD, MAP, and DTCM,  
43 in causing cytotoxicity via both oxidative stress and divalent cation release. These two events are  
44 consistent with cytotoxic mechanisms previously studied for 6HD in models of Parkinsonism in animals.  
45 In these models, 6HD is thought to disrupt mitochondrial Complex I of the electron transport chain,  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 increasing H<sub>2</sub>O<sub>2</sub> production and releasing iron from iron-sulfur centers and heme-iron in Complex I.<sup>28</sup> The  
4 interactions of chelators and of BSO with several neurotoxicants each suggests that disruption of iron  
5 cofactor complexes and depletion of reduced glutathione may be vulnerable toxicity pathways in neurons,  
6 and should be investigated further.  
7  
8  
9  
10

11  
12 A very simple schematic cytotoxicity model is proposed in the Table of Contents graphic to describe 6HD  
13 (or MAP) toxicity *via* both release of iron and by oxidative stress. This schematic shows adaptation of  
14 dLUHMES, by detoxification of 6HD using GSH, as open arrows. 6HD auto-oxidizes to a quinone,  
15 which may be detoxified by conjugation to glutathione.<sup>30</sup> These results do not rule out additional action  
16 of GSH to detoxify peroxides or further oxidation products generated by 6HD disruption of the  
17 mitochondrial electron transport chain, since supplementation with GSH or NAC only partially mitigated  
18 toxicity of 6HD. Modulator BSO inhibits GSH synthesis, and thus blocks detoxification. The solid arrow  
19 in the Table of Contents schematic illustrates 6HD liberation of free Fe<sup>2+</sup> from Complex I of the  
20 mitochondrial electron transport chain, enabling Fe<sup>2+</sup> catalyzation of hydroxyl radical formation which  
21 oxidizes lipids and mitochondrial DNA.<sup>28</sup> Certain chelators, DMPS and PHE, partially protected  
22 dLUHMES cells from 6HD presumably by tying up free Fe<sup>2+</sup>.<sup>28</sup>  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35

## 36 CONCLUSIONS

37  
38  
39 The Toxmatrix method was effective in identifying toxicant:modulator interactions for a variety  
40 of toxicants. This method assayed 9 or 10 concentration levels of each toxicant to identify a  
41 quantitative “tipping point”, the toxicant concentration that overwhelmed the cells ability to  
42 adapt; then examined each modulator for the ability to shift the IC<sub>50</sub> left or right. Each interaction  
43 yielded a hypothesis indicating crosstalk between a toxicity pathway perturbed by the toxicant,  
44 with the pathway addressed by the primary target of the modulator. A variety of  
45 toxicants:modulator interactions identified in these studies represent hypotheses that merit  
46 further investigation, such as the role of sirtuin activator SRT1720 in promoting toxicity of  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

1  
2  
3 MPP+. Given that both activators and inhibitors of sirtuins are proposed as therapeutic  
4  
5 treatments for PD<sup>23</sup>, further interaction studies may facilitate selection of drugs for clinical study.  
6  
7 Several such hypotheses were confirmed for 6HD and MAP toxicity, establishing the liberation  
8  
9 of Fe<sup>2+</sup>, and detoxification by GSH for both toxicants. Modulators such as BSO, DMPS, and  
10  
11 necrostatins were particularly useful because they perturbed toxicity pathways impacting a  
12  
13 variety of toxicants. These results should enable selection of an optimized set of modulators to  
14  
15 improve the efficiency of future Toxmatrix studies by our lab and by others. Thus, this study  
16  
17 harnesses and extends qHTS techniques and analytical tools to demonstrate a new method to  
18  
19 address mechanistic toxicology. Using this new Toxmatrix method, we demonstrate systematic  
20  
21 identification of pathways of toxicity and pathways of cellular adaptation. This is a pathway-  
22  
23 based approach which supports the “Toxicology in the 21<sup>st</sup> century” testing strategy by  
24  
25 identifying molecular changes which alter an adverse outcome. By this approach, chemicals that  
26  
27 are identified as “actives” from quantitative high-throughput screening for stress pathways, can  
28  
29 be further studied to identify toxicity pathways and mechanisms.  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

## ACKNOWLEDGMENT

We are grateful to Craig Thomas, NCATS, for constructive discussion and contribution of thirty-five modulator compounds.

## FUNDING SOURCES

This research was supported by the Tox21 Consortium and the NIH Intramural Research Program of NCATS.

## ABBREVIATIONS

**Toxicants:** 6HD, 6-hydroxydopamine; CAP, Captan; CLM, Chlorambucil; CR2, Sodium dichromate; DGX, Digoxin; DTCM, Dithiocyanatomethane; HCP, Hexachlorophene; MALG, Malachite green oxalate; MAP, 4-(Methylamino)phenol hemisulfate salt; MHG, Methyl mercury (II); MPP, 1-methyl-4-phenylpyridinium; MPTP, 1-methyl-4-phenyl-1,2,3,6-tetrahydropyridine; OTZ, 2-octyl-4-isothiazolin-3-one; TUNC, Tunicamycin; VIN, Vincristine sulfate

**Modulators:** BSO, buthionine sulphoximine; DMPS, 2,3-dimercapto-1-propanesulfonic acid; GSH, reduced glutathione; L-NAME, NG-Nitroarginine Methyl Ester; L-NMMA, NG-Methyl-L-arginine acetate salt; NAC, N-acetyl-L-cysteine; PHE, o-Phenanthroline; and TPEN, (N,N,N',N'-tetrakis(2-pyridylmethyl)-ethylenediamine

## SUPPORTING INFORMATION

A full listing of the 70 modulator chemicals used in this work, SMILE identifiers, and the expected target of each, is provided. This material is freely available via the Internet at <http://pubs.acs.org>.

## REFERENCES

- 1  
2  
3  
4  
5 (1) Xia, M., Huang, R., Witt, K. L., Southall, N., Fostel, J., Cho, M. H., Jadhav, A., Smith, C. S., Inglese, J.,  
6 Portier, C. J., Tice, R. R., and Austin, C. P. (2008) Compound cytotoxicity profiling using quantitative  
7 high-throughput screening. *Environ Health Perspect* 116, 284-291.
- 8 (2) Attene-Ramos, M. S., Miller, N., Huang, R., Michael, S., Itkin, M., Kavlock, R. J., Austin, C. P., Shinn, P.,  
9 Simeonov, A., Tice, R. R., and Xia, M. (2013) The Tox21 robotic platform for the assessment of  
10 environmental chemicals--from vision to reality. *Drug Discov Today* 18, 716-723.
- 11 (3) Huang, R., Xia, M., Sakamuru, S., Zhao, J., Shahane, S. A., Attene-Ramos, M., Zhao, T., Austin, C. P.,  
12 and Simeonov, A. (2016) Modelling the Tox21 10 K chemical profiles for in vivo toxicity prediction  
13 and mechanism characterization. *Nat Commun* 7, 10425.
- 14 (4) Huang, R., Xia, M., Cho, M. H., Sakamuru, S., Shinn, P., Houck, K. A., Dix, D. J., Judson, R. S., Witt, K.  
15 L., Kavlock, R. J., Tice, R. R., and Austin, C. P. (2011) Chemical genomics profiling of environmental  
16 chemical modulation of human nuclear receptors. *Environ Health Perspect* 119, 1142-1148.
- 17 (5) Yilmaz, C., and Ozcengiz, G. (2017) Antibiotics: Pharmacokinetics, toxicity, resistance and multidrug  
18 efflux pumps. *Biochem Pharmacol* 133, 43-62.
- 19 (6) Schildknecht, S., Poltl, D., Nagel, D. M., Matt, F., Scholz, D., Lotharius, J., Schmiegl, N., Salvo-Vargas,  
20 A., and Leist, M. (2009) Requirement of a dopaminergic neuronal phenotype for toxicity of low  
21 concentrations of 1-methyl-4-phenylpyridinium to human cells. *Toxicol Appl Pharmacol* 241, 23-35.
- 22 (7) Bhattacharya, S., Zhang, Q., Carmichael, P. L., Boekelheide, K., and Andersen, M. E. (2011) Toxicity  
23 testing in the 21 century: defining new risk assessment approaches based on perturbation of  
24 intracellular toxicity pathways. *PLoS One* 6, e20887.
- 25 (8) Smirnova, L., Harris, G., Leist, M., and Hartung, T. (2015) Cellular resilience. *Altex* 32, 247-260.
- 26 (9) Tanner, C. M., Goldman, S. M., Ross, G. W., and Grate, S. J. (2014) The disease intersection of  
27 susceptibility and exposure: chemical exposures and neurodegenerative disease risk. *Alzheimers*  
28 *Dement* 10, S213-225.
- 29 (10) Cheung, Y. T., Lau, W. K., Yu, M. S., Lai, C. S., Yeung, S. C., So, K. F., and Chang, R. C. (2009) Effects of  
30 all-trans-retinoic acid on human SH-SY5Y neuroblastoma as in vitro model in neurotoxicity research.  
31 *Neurotoxicology* 30, 127-135.
- 32 (11) Lotharius, J., Falsig, J., van Beek, J., Payne, S., Dringen, R., Brundin, P., and Leist, M. (2005)  
33 Progressive degeneration of human mesencephalic neuron-derived cells triggered by dopamine-  
34 dependent oxidative stress is dependent on the mixed-lineage kinase pathway. *J Neurosci* 25, 6329-  
35 6342.
- 36 (12) Tong, Z. B., Hogberg, H., Kuo, D., Sakamuru, S., Xia, M., Smirnova, L., Hartung, T., and Gerhold, D.  
37 (2017) Characterization of three human cell line models for high-throughput neuronal cytotoxicity  
38 screening. *J Appl Toxicol* 37, 167-180.
- 39 (13) Shin, S., and Vemuri, M. (2010) Culture and Differentiation of Human Neural Stem Cells. *Springer*  
40 *Protoc Hand*, 51-73.
- 41 (14) Wolpaw, A. J., Shimada, K., Skouta, R., Welsch, M. E., Akavia, U. D., Pe'er, D., Shaik, F., Bulinski, J. C.,  
42 and Stockwell, B. R. (2011) Modulatory profiling identifies mechanisms of small molecule-induced  
43 cell death. *Proc Natl Acad Sci U S A* 108, E771-780.
- 44 (15) Xia, M., Huang, R., Witt, K. L., Southall, N., Fostel, J., Cho, M. H., Jadhav, A., Smith, C. S., Inglese, J.,  
45 Portier, C. J., Tice, R. R., and Austin, C. P. (2008) Compound cytotoxicity profiling using quantitative  
46 high-throughput screening. *Environ Health Perspect* 116, 284-291.
- 47 (16) Wang, Y., and Huang, R. (2016) Correction of Microplate Data from High Throughput Screening, In  
48 *High-Throughput Screening Assays in Toxicology* (Zhu, H., and Xia, M., Eds.), Humana Press.
- 49 (17) Hill, A. V. (1910) The possible effects of the aggregation of the molecules of haemoglobin on its  
50 dissociation curves. *J. Physiol. (London)* 40, 4-7.
- 51  
52  
53  
54  
55  
56  
57  
58  
59  
60

- 1  
2  
3 (18)Wang, Y., Jadhav, A., Southal, N., Huang, R., and Nguyen, D. T. (2010) A grid algorithm for high  
4 throughput fitting of dose-response curve data. *Curr Chem Genomics* 4, 57-66.
- 5 (19)Inglese, J., Auld, D. S., Jadhav, A., Johnson, R. L., Simeonov, A., Yasgar, A., Zheng, W., and Austin, C.  
6 P. (2006) Quantitative high-throughput screening: a titration-based approach that efficiently  
7 identifies biological activities in large chemical libraries. *Proc Natl Acad Sci U S A* 103, 11473-11478.
- 8 (20)Zheng, W., Degtarev, A., Hsu, E., Yuan, J., and Yuan, C. (2008) Structure-activity relationship study of  
9 a novel necroptosis inhibitor, necrostatin-7. *Bioorg Med Chem Lett* 18, 4932-4935.
- 10 (21)O'Donnell, M. A., and Ting, A. T. (2010) Chronicles of a death foretold: dual sequential cell death  
11 checkpoints in TNF signaling. *Cell Cycle* 9, 1065-1071.
- 12 (22)Boll, B., Eltaib, F., Reiners, K. S., von Tresckow, B., Tawadros, S., Simhadri, V. R., Burrows, F. J.,  
13 Lundgren, K., Hansen, H. P., Engert, A., and von Strandmann, E. P. (2009) Heat shock protein 90  
14 inhibitor BIIB021 (CNF2024) depletes NF-kappaB and sensitizes Hodgkin's lymphoma cells for natural  
15 killer cell-mediated cytotoxicity. *Clin Cancer Res* 15, 5108-5116.
- 16 (23)Min, S. W., Sohn, P. D., Cho, S. H., Swanson, R. A., and Gan, L. (2013) Sirtuins in neurodegenerative  
17 diseases: an update on potential mechanisms. *Front Aging Neurosci* 5, 53.
- 18 (24)Rowlands, B. D., Lau, C. L., Ryall, J. G., Thomas, D. S., Klugmann, M., Beart, P. M., and Rae, C. D.  
19 (2015) Silent information regulator 1 modulator resveratrol increases brain lactate production and  
20 inhibits mitochondrial metabolism, whereas SRT1720 increases oxidative metabolism. *J Neurosci Res*  
21 93, 1147-1156.
- 22 (25)Sabbatini, P., Rowand, J. L., Groy, A., Korenchuk, S., Liu, Q., Atkins, C., Dumble, M., Yang, J.,  
23 Anderson, K., Wilson, B. J., Emmitte, K. A., Rabindran, S. K., and Kumar, R. (2009) Antitumor activity  
24 of GSK1904529A, a small-molecule inhibitor of the insulin-like growth factor-I receptor tyrosine  
25 kinase. *Clin Cancer Res* 15, 3058-3067.
- 26 (26)Renault, T. T., Elkholi, R., Bharti, A., and Chipuk, J. E. (2014) B cell lymphoma-2 (BCL-2) homology  
27 domain 3 (BH3) mimetics demonstrate differential activities dependent upon the functional  
28 repertoire of pro- and anti-apoptotic BCL-2 family proteins. *J Biol Chem* 289, 26481-26491.
- 29 (27)Shimizu, E., Hashimoto, K., Komatsu, N., and Iyo, M. (2002) Roles of endogenous glutathione levels  
30 on 6-hydroxydopamine-induced apoptotic neuronal cell death in human neuroblastoma SK-N-SH  
31 cells. *Neuropharmacology* 43, 434-443.
- 32 (28)Weinreb, O., Mandel, S., Youdim, M. B., and Amit, T. (2013) Targeting dysregulation of brain iron  
33 homeostasis in Parkinson's disease by iron chelators. *Free Radic Biol Med* 62, 52-64.
- 34 (29)Drew, R., and Miners, J. O. (1984) The effects of buthionine sulphoximine (BSO) on glutathione  
35 depletion and xenobiotic biotransformation. *Biochem Pharmacol* 33, 2989-2994.
- 36 (30)Izumi, Y., Sawada, H., Sakka, N., Yamamoto, N., Kume, T., Katsuki, H., Shimohama, S., and Akaike, A.  
37 (2005) p-Quinone mediates 6-hydroxydopamine-induced dopaminergic neuronal death and ferrous  
38 iron accelerates the conversion of p-quinone into melanin extracellularly. *J Neurosci Res* 79, 849-  
39 860.
- 40 (31)Tirmenstein, M. A., Hu, C. X., Scicchitano, M. S., Narayanan, P. K., McFarland, D. C., Thomas, H. C.,  
41 and Schwartz, L. W. (2005) Effects of 6-hydroxydopamine on mitochondrial function and glutathione  
42 status in SH-SY5Y human neuroblastoma cells. *Toxicol In Vitro* 19, 471-479.
- 43 (32)Aposhian, H. V., Maiorino, R. M., Rivera, M., Bruce, D. C., Dart, R. C., Hurlbut, K. M., Levine, D. J.,  
44 Zheng, W., Fernando, Q., Carter, D., and et al. (1992) Human studies with the chelating agents,  
45 DMPS and DMSA. *J Toxicol Clin Toxicol* 30, 505-528.
- 46 (33)Carranza-Rosales, P., Guzman-Delgado, N. E., Cruz-Vega, D. E., Balderas-Renteria, I., and Gandolfi, A.  
47 J. (2007) DMPS reverts morphologic and mitochondrial damage in OK cells exposed to toxic  
48 concentrations of HgCl<sub>2</sub>. *Cell Biol Toxicol* 23, 163-176.
- 49 (34)Youdim, M. B. (2013) Multi target neuroprotective and neurorestorative anti-Parkinson and anti-  
50 Alzheimer drugs ladostigil and m30 derived from rasagiline. *Exp Neurobiol* 22, 1-10.
- 51  
52  
53  
54  
55  
56  
57  
58  
59  
60

- 1  
2  
3 (35) Kobayashi, H., Oikawa, S., Umemura, S., Hirose, I., and Kawanishi, S. (2008) Mechanism of metal-  
4 mediated DNA damage and apoptosis induced by 6-hydroxydopamine in neuroblastoma SH-SY5Y  
5 cells. *Free Radic Res* 42, 651-660.  
6  
7 (36) Fine, J. M., Forsberg, A. C., Renner, D. B., Faltsek, K. A., Mohan, K. G., Wong, J. C., Arneson, L. C.,  
8 Crow, J. M., Frey, W. H., 2nd, and Hanson, L. R. (2014) Intranasally-administered deferoxamine  
9 mitigates toxicity of 6-OHDA in a rat model of Parkinson's disease. *Brain Res* 1574, 96-104.  
10 (37) Devos, D., Moreau, C., Devedjian, J. C., Kluz, J., Petrault, M., Laloux, C., Jonneaux, A., Ryckewaert,  
11 G., Garcon, G., Rouaix, N., Duhamel, A., Jissendi, P., Dujardin, K., Auger, F., Ravasi, L., Hopes, L.,  
12 Grolez, G., Firdaus, W., Sablonniere, B., Strubi-Vuillaume, I., Zahr, N., Destee, A., Corvol, J. C., Poltl,  
13 D., Leist, M., Rose, C., Defebvre, L., Marchetti, P., Cabantchik, Z. I., and Bordet, R. (2014) Targeting  
14 chelatable iron as a therapeutic modality in Parkinson's disease. *Antioxid Redox Signal* 21, 195-210.  
15 (38) Martin-Bastida, A., Ward, R. J., Newbould, R., Piccini, P., Sharp, D., Kabba, C., Patel, M. C., Spino, M.,  
16 Connelly, J., Tricta, F., Crichton, R. R., and Dexter, D. T. (2017) Brain iron chelation by deferiprone in  
17 a phase 2 randomised double-blinded placebo controlled clinical trial in Parkinson's disease. *Sci Rep*  
18 7, 1398.  
19  
20  
21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

## Figure Legends

**Table 1.** Modulators, CAS # identifiers, and expected targets for each. Abbreviations: BSO, buthionine sulphoximine; DMPS, 2,3-dimercapto-1-propanesulfonic acid; GSH, reduced glutathione; L-NAME, NG-Nitroarginine Methyl Ester; L-NMMA, NG-Methyl-L-arginine acetate salt; NAC, N-acetyl-L-cysteine; PHE, o-Phenanthroline; and TPEN, (N,N,N', N'-tetrakis(2-pyridylmethyl)-ethylenediamine; IGF1R, Insulin-like Growth Factor 1 Receptor; IR, Insulin Receptor; NOS, Nitric Oxide Synthase (-1, -2, and -3); A full list of the 70 modulators used in screenings is available in Supplementary Table S1.

**Figure 1.** Example plots illustrating interactions between toxicant MAP and indicated modulators in SH-SY5Y cells: sensitization to cytotoxicity (**A**), no interaction (**B**), or protection (**C**). Green spots show that the modulators were not toxic at the three concentrations indicated, except the highest concentration of  $\text{CoCl}_2$  which killed the cells. Black lines indicate toxicant MAP alone concentration-response; Dark blue-, Light blue-, and grey-, lines indicate MAP plus a modulator at low, medium, and high concentrations, respectively.

**Figure 2.** Heatmap showing interactions between toxicants<sup>12</sup> and modulators (columns) that alter cytotoxicity in SH-SY5Y cells. Maximal fold-change values are shown comparing  $\text{IC}_{50}$  value with modulator divided by  $\text{IC}_{50}$  value without modulator. Cell sensitization to a toxicant by a modulator is indicated by a positive  $\text{IC}_{50}$  shift value if it exceeded 3-fold increase or decrease at the highest modulator concentration that did not cause  $\geq 10\%$  cytotoxicity. Red fill indicates sensitization to a toxicant, with increasing color saturation from 2 to 16-fold. Blue fill indicates a negative  $\text{IC}_{50}$  shift hence protection from a toxicant, with increasing color saturation from -2 to -16-fold. White fill indicates no interaction. 32 toxicants were tested in 10 concentrations, or with DMSO vehicle, and with 3 concentrations of each of 70 modulators. Toxicant/modulator combinations that yielded significant interactions, are colored red for increased sensitivity, blue for protection by modulator, or white for no significant interaction. The maximum fold-change shift in  $\text{IC}_{50}$  observed is shown in each colored square.


1  
2  
3 **Figure 3.** Heatmap illustrating toxicant:modulator interactions in three differentiated neuronal models.  
4  
5 Differentiated human cells are labeled “dSH-SY5Y”, “dLUHMES”, or “dNSC”. Among  
6  
7 toxicant:modulator pairs tested, 16 did not show IC<sub>50</sub> shifts exceeding 3-fold for any of the three  
8  
9 differentiated cell types, and are omitted from this Figure. Cell colors are as described in the Figure 2  
10  
11 legend except extended to show 2- to 16-fold IC<sub>50</sub> shifts. Modulator abbreviations and drug targets are  
12  
13 listed in Table 1. “NT” indicates the toxicant did not kill the cells at the highest concentration tested, 100  
14  
15 μM; and “TOX” indicates the toxicant killed the cells at the lowest concentration tested, 15 nM  
16  
17

18 **Figure 4.** Sensitization of dLUHMES cells to 10 μM 6HD (A), or 4 μM MAP (B) by glutathione  
19  
20 synthesis inhibitor, BSO, is partially relieved by 100 μM GSH or 100 μM NAC. BSO alone (open circles)  
21  
22 was moderately cytotoxic to dLUHMES cells; whereas A) 10 μM 6HD + BSO (filled squares) were  
23  
24 highly toxic. Addition of 100 μM NAC (grey filled circles) or 100 μM GSH (open triangles) to 6HD +  
25  
26 BSO or MAP + BSO moderated cytotoxicity.  
27  
28  
29


30 **Figure 5.** Protection of dLUHMES cells from 6HD MAP by PHE, but not by M30. PHE alone (open  
31  
32 circles) was not toxic to dLUHMES cells; whereas M30 alone (open squares) was moderately cytotoxic.  
33  
34 Addition of PHE (filled circles) to 10 μM 6HD (A), or 4 μM MAP (B) moderated cytotoxicity; whereas  
35  
36 addition of M30 (filled squares) increased cytotoxicity.  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

Modulator	CAS #	Target
BIIB021	848695-25-0	HSP90 inhibitor
BSO	5072-26-4	Glutamate Cysteine Ligase inhibitor
CoCl <sub>2</sub>	7646-79-9	HIF1a/hypoxia inducer
DMPS	74-61-3	divalent cation chelator
Eliprodil	119431-25-3	NMDA Antagonist
GSK-1904529A	1089283-49-7	IGF-1R/IR inhibitor
KU0063794	938440-64-3	mTOR inhibitor
L-NAME	50903-99-6	NOS inhibitor
NAC	616-91-1	ROS scavenger
necrostatin 1	4311-88-0	RIP1K inhibitor
necrostatin 5	337349-54-9	RIP1K inhibitor
necrostatin 7	351062-08-3	necroptosis inhibitor
Nicotinamide	98-92-0	SIRT1/2/3/6 inhibitor
o-Phenanthroline	3248-05-3	Chelator of divalent cations
Pazopanib	444731-52-6	VEGFR1/2/3 inhibitor
Saquinavir	149845-06-7	HIV protease- and CYP3A4 inhibitor
SP-600125	129-56-6	JNK1, JNK2, JNK3 inhibitor
SRT1720	1001645-58-4	SIRT1 activator
TPEN	16858-02-9	divalent cation chelator
TW-37	877877-35-5	BCL2-binding inhibitor
WYE-354	1062169-56-5	mTOR inhibitor

**Table 1.** Modulators, CAS # identifiers, and expected targets for each. Abbreviations: BSO, buthionine sulphoximine; DMPS, 2,3-dimercapto-1-propanesulfonic acid; GSH, reduced glutathione; L-NAME, NG-Nitroarginine Methyl Ester; L-NMMA, NG-Methyl-L-arginine acetate salt; NAC, N-acetyl-L-cysteine; PHE, o-Phenanthroline; and TPEN, (N,N,N', N'-tetrakis(2-pyridylmethyl)-ethylenediamine; IGF1R, Insulin-like Growth Factor 1 Receptor; IR, Insulin Receptor; NOS, Nitric Oxide Synthase (-1, -2, and -3); A full list of the 70 modulators used in screenings is available in Supplementary Table S1.





**Figure 1.** Example plots illustrating interactions between toxicant MAP and indicated modulators in SH-SY5Y cells: sensitization to cytotoxicity (a), no interaction (b), or protection (c). Green spots show that the modulators were not toxic at the three concentrations indicated, except the highest concentration of  $\text{CoCl}_2$  which killed the cells. Black lines indicate toxicant MAP alone concentration-response; Dark blue-, Light blue-, and grey-, lines indicate MAP plus a modulator at low, medium, and high concentrations, respectively. In order to illustrate variance the 12 data points are shown for the toxicant-only samples from the same microplate, and the 8 data points are shown for the modulator-only controls from the same microplate.





**Figure 2.** Heatmap showing interactions between toxicants<sup>12</sup> and modulators (columns) that alter cytotoxicity in SH-SY5Y cells. Maximal fold-change values are shown comparing IC<sub>50</sub> value with modulator divided by IC<sub>50</sub> value without modulator. Cell sensitization to a toxicant by a modulator is indicated by a positive IC<sub>50</sub> shift value if it exceeded 3-fold increase or decrease at the highest modulator concentration that did not cause  $\geq 10\%$  cytotoxicity. Red fill indicates sensitization to a toxicant, with increasing color saturation from 2- to 16-fold. Blue fill indicates a negative IC<sub>50</sub> shift hence protection from a toxicant, with increasing color saturation from negative 2- to 16-fold. White fill indicates no interaction. 32 toxicants were tested in 10 concentrations, or with DMSO vehicle, and with 3 concentrations of each of 70 modulators. Toxicant/modulator combinations that yielded significant interactions, are colored red for increased sensitivity, blue for protection by modulator, or white for no significant interaction. The maximum fold-change shift in IC<sub>50</sub> observed is shown in each colored square.

Toxicant	Modulator	dSH-SY5Y		dLUHMES		dNSC	
		24h	48h	24h	48h	24h	48h
MAP	CoCl <sub>2</sub>			4.5	13		6.3
MAP	DMPS	-3.2	-2.8	-2.5	-2.4	-3.1	
MAP	N-acetyl-cysteine	-3.8	-3.9				
6HD	BSO		4.3	910			72
6HD	DMPS	-2.9		-6.1	-3.5	-3.8	
6HD	Necrostatin 7		25	74	23		
6HD	WYE-354					4.9	21
CAP	DMPS	-11	-8.7	-8.3		-4.2	-3.3
CLM	CNF-2024		3.2	4.5	3.6	2.4	2.5
DTCM	DMPS	-18	-6.4	-10	-8.6	-2.6	
DTCM	N-acetyl-cysteine	-9.7		-9.4	-6.2		-2.3
MPP	SRT1720	4.8		4.2	3.0		3.8
MALG	GSK-1904529A	3.5					
MALG	Necrostatin 7		3.4				4.3
MHG	DMPS	-2.7		-2.5	-9.0	-7.3	
MHG	N-acetyl-cysteine	-2.1	-3.0	-2.2	-2.5		-2.3
MHG	Necrostatin 7		10	4.6		2.2	5.9
MHG	Nicotinamide				-4.0		
TUNC	TW-37	10				NT	
VIN	SP-600125			-4.4		-11	TOX

**Figure 3.** Heatmap illustrating toxicant:modulator interactions in three differentiated neuronal models. Differentiated human cells are labeled “dSH-SY5Y”, “dLUHMES”, or “dNSC”. Among toxicant:modulator pairs tested, 16 did not show IC<sub>50</sub> shifts exceeding 3-fold for any of the three differentiated cell types, and are omitted from this Figure. Cell colors are as described in the Figure 2 legend except extended to show 2- to 16-fold IC<sub>50</sub> shifts. Modulator abbreviations and drug targets are listed in Table 1. “NT” indicates the toxicant did not kill the cells at the highest concentration tested, 100 μM; and “TOX” indicates the toxicant killed the cells at the lowest concentration tested, 15 nM.



**Figure 4.** Sensitization of dLUHMES cells to 10  $\mu\text{M}$  6HD (a), or 4  $\mu\text{M}$  MAP (b) by glutathione synthesis inhibitor, BSO, is partially relieved by 100  $\mu\text{M}$  GSH or 100  $\mu\text{M}$  NAC. BSO alone (open circles) was moderately cytotoxic to dLUHMES cells; whereas A) 10  $\mu\text{M}$  6HD + BSO (filled squares) were highly toxic. Addition of 100  $\mu\text{M}$  NAC (grey filled circles) or 100  $\mu\text{M}$  GSH (open triangles) to 6HD + BSO or MAP + BSO moderated cytotoxicity.



50  
51  
52  
53  
54  
55  
56  
57  
58  
59  
60

**Figure 5.** Protection of dLUHMES cells from 6HD MAP by PHE, but not by M30. PHE alone (open circles) was not toxic to dLUHMES cells; whereas M30 alone (open squares) was moderately cytotoxic. Addition of PHE (filled circles) to 10  $\mu$ M 6HD (A), or 4  $\mu$ M MAP (B) moderated cytotoxicity; whereas addition of M30 (filled squares) increased cytotoxicity.